Lab 3: Monod cell growth kinetics in batch reactor and chemostat

A microorganism is known to follow Monod kinetics with the following parameters;

$$\mu_{\text{max}} = 1 \text{ h}^{-1}$$
; $Y_{xs} = 0.5$; and $K_s = 0.1 \text{ g/l}$

Batch culture

- a. If the microorganism is grown in batch culture, obtain a plot for the change in cell concentration (X₀ = 0.1 g/l; S₀ = 5 g/l) with time by employing numerical expression using Monod's equation and cell growth equation as discussed in class.
- b. Use the analytical solution for the problem discussed in part a and obtain a plot of change in cell concentration with time. Compare the results of part a and b.

Chemostat

c. If the same microorganism is grown in a chemostat with the following (X₀ = 0.1 g/l; S₀ = 5 g/l), obtain plots for the change in cell and substrate concentration with time till steady state at the following dilution rates (D)

$$[0.2 h^{-1}; 0.4 h^{-1}; 0.6 h^{-1}; 0.8 h^{-1}$$
 and $1 h^{-1})$.

- d. Based on the steady state values calculated in part b, plot the following at the 5 dilution rates D.
 - Substrate concentration in chemostat, S
 - ii. Steady state cell concentration, X
 - Output of cells from chemostat, DX